




1



2




4



5



6



7

Recap of Day 1

- Data Quality
- Acceleration metrics or features
- Adapt GGIR to different study protocols
- Descriptive measures and physical activity

8

Configuration Parameters Vignette



9

Day Segment Analysis Vignette



10

Assignment 1

1. Open RStudio and an empty script
2. Create a GGIR function call
3. Define datadir and outputdir
 - Tip 1: datadir should specify the path to out demo file
 - Tip 2: outputdir should be an existing folder (different to datadir)
4. Define mode to run GGIR parts 1 and 2
5. Make sure you only analyse data from the first midnight onwards
6. We are only interested in the analysis of the first 3 days.
7. Run the script via the source button
8. Advanced: Look up the output and meaning of variables
9. Optional: Try to run GGIR parts 1 and 2 on your own data


11

Part 2 Output

www.accelting.com

12

The GGIR pipeline



13

Part 2

Output

- meta
 - csv
 - IDo1.csv
 - ...
- results
 - QC
 - data_quality_report.csv
 - plot_to_check_data_quality.pdf
 - part2_daysummary.csv
 - part2_daysummary_longformat.csv
 - part2_summary.csv



14

Part 2

Output

- meta
 - csv
 - IDo1.csv
 - ...
- results
 - QC
 - data_quality_report.csv
 - plot_to_check_data_quality.pdf
 - part2_daysummary.csv**
 - part2_daysummary_longformat.csv**
 - part2_summary.csv



15

Output from Part 2

Day-level features (wide)

do.report = 2

ID	filename	calendar_date	bodylocation	N valid hours	N hours	weekday	Measurement day	qwindow_timestamps	qwindow_name
11	011_45400.cwa	2019-03-28T00:00:00+0100	not extracted	24	24	Wednesday	2	0_8_24	0_8_24
11	011_45400.cwa	2019-03-29T00:00:00+0100	not extracted	24	24	Thursday	3	0_8_24	0_8_24
11	011_45400.cwa	2019-03-30T00:00:00+0100	not extracted	24	24	Friday	4	0_8_24	0_8_24
11	011_45400.cwa	2019-03-31T00:00:00+0100	not extracted	23	23	Saturday	5	0_8_24	0_8_24

16

Output from Part 2

Day-level features (wide)

do.report = 2

part2_daysummary.csv

ID	Measurement day	qwindow_timestamps	mean_ENMO_mg_0-24hr	MVPA_ESS_T100_ENMO_0-24hr	mean_ENMO_mg_0-8hr	MVPA_ESS_T100_ENMO_0-8hr	mean_ENMO_mg_8-24hr	MVPA_ESS_T100_ENMO_8-24hr
11	2	0_8_24	10.291	146.833	5.368	2.107	17.762	144.667
11	3	0_8_24	10.099	151.833	4.026	3.417	22.136	148.417
11	4	0_8_24	10.233	170.417	7.903	8.5	13.386	161.917
11	5	0_8_24	11.085	41.25	7.383	4.167	18.932	37.083

17

Output from Part 2

Day-level features (long)

do.report = 2 # Only if different qwindows are defined (e.g., qwindow = c(0, 8, 24))

part2_daysummary_longformat.csv

ID	filename	calendar_date	bodylocation	N_valid_hours	N_hours	N_valid_hours_in_window	N_hours_in_window	weekday	Measurement day	qwindow_timestamps	qwindow_name
11	011_45400.cwa	2019-03-28T00:00:00+0100	not extracted	24	24	8	8	Wednesday	2	00:00:00:00	0-8hr
11	011_45400.cwa	2019-03-29T00:00:00+0100	not extracted	24	24	8	8	Thursday	3	00:00:00:00	0-8hr
11	011_45400.cwa	2019-03-30T00:00:00+0100	not extracted	24	24	8	8	Friday	4	00:00:00:00	0-8hr
11	011_45400.cwa	2019-03-31T00:00:00+0100	not extracted	23	23	8	8	Saturday	5	00:00:00:00	0-8hr
11	011_45400.cwa	2019-03-28T00:00:00+0100	not extracted	24	24	16	16	Wednesday	2	08:00:00:00	8-24hr
11	011_45400.cwa	2019-03-29T00:00:00+0100	not extracted	24	24	16	16	Thursday	3	08:00:00:00	8-24hr
11	011_45400.cwa	2019-03-30T00:00:00+0100	not extracted	24	24	16	16	Friday	4	08:00:00:00	8-24hr
11	011_45400.cwa	2019-03-31T00:00:00+0100	not extracted	23	23	15	15	Saturday	5	08:00:00:00	8-24hr
11	011_45400.cwa	2019-03-28T00:00:00+0100	not extracted	24	24			Wednesday	2	00:00:00:00	0-24hr
11	011_45400.cwa	2019-03-29T00:00:00+0100	not extracted	24	24			Thursday	3	00:00:00:00	0-24hr
11	011_45400.cwa	2019-03-30T00:00:00+0100	not extracted	24	24			Friday	4	00:00:00:00	0-24hr
11	011_45400.cwa	2019-03-31T00:00:00+0100	not extracted	23	23			Saturday	5	00:00:00:00	0-24hr

18

Output from Part 2

Day-level features (long)

do.report = 2 # Only if different qwindows are defined (e.g., qwindow = c(0, 8, 24))

part2_daysummary_longformat.csv

ID	filename	calendar_date	bodylocation	N_valid_hours	N_hours	N_valid_hours_in_window	N_hours_in_window	weekday	Measurement day	qwindow_timestamps	qwindow_name
11	011_45400.cwa	2019-03-28T00:00:00+0100	not extracted	24	24	8	8	Wednesday	2	00:00:00:00	0-8hr
11	011_45400.cwa	2019-03-29T00:00:00+0100	not extracted	24	24	8	8	Thursday	3	00:00:00:00	0-8hr
11	011_45400.cwa	2019-03-30T00:00:00+0100	not extracted	24	24	8	8	Friday	4	00:00:00:00	0-8hr
11	011_45400.cwa	2019-03-31T00:00:00+0100	not extracted	23	23	8	8	Saturday	5	00:00:00:00	0-8hr
11	011_45400.cwa	2019-03-28T00:00:00+0100	not extracted	24	24	16	16	Wednesday	2	08:00:00:00	8-24hr
11	011_45400.cwa	2019-03-29T00:00:00+0100	not extracted	24	24	16	16	Thursday	3	08:00:00:00	8-24hr
11	011_45400.cwa	2019-03-30T00:00:00+0100	not extracted	24	24	16	16	Friday	4	08:00:00:00	8-24hr
11	011_45400.cwa	2019-03-31T00:00:00+0100	not extracted	23	23	15	15	Saturday	5	08:00:00:00	8-24hr
11	011_45400.cwa	2019-03-28T00:00:00+0100	not extracted	24	24			Wednesday	2	00:00:00:00	0-24hr
11	011_45400.cwa	2019-03-29T00:00:00+0100	not extracted	24	24			Thursday	3	00:00:00:00	0-24hr
11	011_45400.cwa	2019-03-30T00:00:00+0100	not extracted	24	24			Friday	4	00:00:00:00	0-24hr
11	011_45400.cwa	2019-03-31T00:00:00+0100	not extracted	23	23			Saturday	5	00:00:00:00	0-24hr

19

Output from Part 2

Day-level features (long)

do_report = 2 # default = c(2, 4, 5) - Only if different qwindows are defined (e.g., qwindow = c(0, 8, 24))

part2_daysummary_longformat.csv

ID	qwindow_timestamps	qwindow_name	mean_ENMO_mg	MVPA_ESS_T100_ENMO	MVPA_E1M_T100_ENMO	MVPA_ESM_T100_ENMO	MVPA_ESS_B1MB0%_T100_ENMO	MVPA_ESS_B5MB0%_T100_ENMO	MVPA_ESS_B10MB0%_T100_ENMO
11	00:00-8:00	0-8hr	5.388	2.167	0	0	0	0	0
11	00:00-8:00	0-8hr	4.024	2.437	0	0	2.167	0	0
11	00:00-8:00	0-8hr	7.393	8.5	0	0	2.917	0	0
11	00:00-8:00	0-8hr	7.393	4.167	4	0	5.5	0	0
11	08:00-24:00	8-24hr	22.762	54.667	155	140	103.917	54.25	27.417
11	08:00-24:00	8-24hr	22.136	48.417	37	10	15.083	0	0
11	08:00-24:00	8-24hr	53.396	161.917	272	160	114.167	103.833	96.833
11	08:00-24:00	8-24hr	18.932	37.083	19	0	11.417	5.167	0
11	00:00-24:00	0-24hr	10.297	146.833	155	140	103.917	54.25	27.417
11	00:00-24:00	0-24hr	16.099	51.833	40	10	16.25	0	0
11	00:00-24:00	0-24hr	38.232	170.417	181	170	117.083	103.833	96.833
11	00:00-24:00	0-24hr	15.083	81.25	33	0	12.917	5.167	0

20

Part 2

Output

- meta
 - csv
 - ID01.csv
 - ...
- results
 - QC
 - data_quality_report.csv
 - plot_to_check_data_quality.pdf
 - part2_daysummary.csv
 - part2_daysummary_longformat.csv
 - part2_summary.csv



21

Output from Part 2

Person-level features

do_report = 2

part2_summary.csv

ID	device_sn	bodylocation	filename	start_time	startdate	samplefreq	device	clipping_score	meas_dur_days	complete_24hrcycle
11	45400	not extracted	011_45400.cwa	2019-03-27T00:00:00+0100	Tuesday	50	activity	0	6.99	1

All days

ID	N valid Weddays	N valid Wkdays	AD_mean_ENMO_mg_0-24hr	AD_MVPA_ESS_T100_ENMO_0-24hr	AD_mean_ENMO_mg_0-8hr	AD_MVPA_ESS_T100_ENMO_0-8hr	AD_mean_ENMO_mg_8-24hr	AD_MVPA_ESS_T100_ENMO_8-24hr
11	2	5	28.809	94.81	6.222	4.083	40.102	90.726

Weekdays

ID	N valid Weddays	N valid Wkdays	WD_mean_ENMO_mg_0-24hr	WD_MVPA_ESS_T100_ENMO_0-24hr	WD_mean_ENMO_mg_0-8hr	WD_MVPA_ESS_T100_ENMO_0-8hr	WD_mean_ENMO_mg_8-24hr	WD_MVPA_ESS_T100_ENMO_8-24hr
11	2	5	33.197	111.05	5.471	4	47.06	107.05

Weekend days

ID	N valid Weddays	N valid Wkdays	WE_mean_ENMO_mg_0-24hr	WE_MVPA_ESS_T100_ENMO_0-24hr	WE_mean_ENMO_mg_0-8hr	WE_MVPA_ESS_T100_ENMO_0-8hr	WE_mean_ENMO_mg_8-24hr	WE_MVPA_ESS_T100_ENMO_8-24hr
11	2	5	17.838	54.208	8.099	4.292	22.707	48.917

22

Part 2

Anything else?

- Descriptives of behaviour
 - Intensity gradient (Rowlands et al)
 - MX metrics
 - M5 & L5 metrics
 - Interdaily stability (IS)
 - Intradaily variability (IV)
- Output
 - Weighted weekdays and weekend days



23

Part 2

Output

- meta
 - csv
 - ID01.csv
 - ...
- results
 - QC
 - data_quality_report.csv
 - plot_to_check_data_quality.pdf
 - part2_daysummary.csv
 - part2_daysummary_longformat.csv
 - part2_summary.csv



24

Output from Part 2

Data quality reports

output_activity > results > QC		
Name	Status	Date modified
data_quality_report.csv	✓	6/2/2022 9:47 AM
plots_to_check_data_quality_1.pdf	✓	6/2/2022 9:47 AM


Filename	cal.error.start	cal.error.end	QCmessage
meta_011_45400.cwa.RData	0.12640	0.00466	recalibration done, no problems detected
meta_013_42151.cwa.RData	0.05089	0.00286	recalibration done, no problems detected
meta_015_44944.cwa.RData	0.08323	0.00355	recalibration done, no problems detected
meta_017_42151.cwa.RData	0.03493	0.00270	recalibration done, no problems detected

25

Output from Part 2


Data quality reports

do part2pdf = TRUE # default = TRUE




26

The GGIR pipeline (output folder)



27

The GGIR pipeline (output folder)



28

Part 2

Output

- meta
 - CSV
 - IDO1.csv
 - ...
- results
 - QC
 - data_quality_report.csv
 - plot_to_check_data_quality.pdf
 - part2_daysummary.csv
 - part2_daysummary_longformat.csv
 - part2_summary.csv



29

Output from Part 2

Epoch-level features

epochvalues2csv = TRUE # default = FALSE!

Name	Status	Date modified
011_45400_cwa.cwa.RData.csv	Ⓞ	6/2/2022 9:46 AM
011_42151_cwa.cwa.RData.csv	Ⓞ	6/2/2022 9:46 AM
011_44944_cwa.cwa.RData.csv	Ⓞ	6/2/2022 9:46 AM
011_42151_cwa.cwa.RData.csv	Ⓞ	6/2/2022 9:46 AM

Signal metrics only!

For time series of sleep and bout classification, see GGIR part 5 (as discussed on day 3)

	A	B	C
1	timestamp	anglez	ENMO
2	2019-03-27T00:00:00+0100	-25.347	0.0022
3	2019-03-27T00:00:05+0100	-25.3816	0.0017
4	2019-03-27T00:00:10+0100	-25.342	0.0025
5	2019-03-27T00:00:15+0100	-24.8968	0.0027
6	2019-03-27T00:00:20+0100	-25.3328	0.0027
7	2019-03-27T00:00:25+0100	-25.4124	0.0024
8	2019-03-27T00:00:30+0100	-25.2169	0.002
9	2019-03-27T00:00:35+0100	-25.7494	0.0023
10	2019-03-27T00:00:40+0100	-25.5495	0.0024
11	2019-03-27T00:00:45+0100	-26.3225	0.0021
12	2019-03-27T00:00:50+0100	-25.4574	0.0022
13	2019-03-27T00:00:55+0100	-25.2172	0.0028
14	2019-03-27T00:01:00+0100	-24.6724	0.0025
15	2019-03-27T00:01:05+0100	-25.0957	0.0023
16	2019-03-27T00:01:10+0100	-24.5152	0.0026
17	2019-03-27T00:01:15+0100	-24.996	0.0033
18	2019-03-27T00:01:20+0100	-26.0688	0.0019
19	2019-03-27T00:01:25+0100	-25.278	0.0028
20	2019-03-27T00:01:30+0100	-25.5993	0.0023
21	2019-03-27T00:01:35+0100	-25.8046	0.0024

30

Assignment 2 (build on assignment 1 script)

1. Derive MX metrics (Rowlands et al.): M120, M60, M30, M5
2. Derive intensity levels: 0, 50, 100, 200, 400, 800
3. Derive the intensity gradient
4. Derive the MVPA thresholds: 100 and 140
5. Advanced: Calculate these variables over the following Windows:
 - From 0 am to 7am
 - From 7am to 10pm
 - From 10pm to 0am next day

31



32

Polysomnography

Lab-based

- Reference method for sleep assessment
- Combination of monitors:
 - EEG
 - Heart rate
 - Gases exchange
 - Blood oxygen levels
 - Others
- Specialist go over the signal and classify 30-sec epochs into sleep stages

33

Accelerometer-based sleep assessment

Free living

- Challenge 1:** distinguish sleep, wake, and non-wear
 - Lack of movement
 - Lack of postural change

- Term: Sustained Inactivity Bout, abbreviated as SIB

34

Accelerometer-based sleep assessment

Free living

- Challenge 2:** separate daytime and nighttime (**Sleep Period Time**)
 - Automatic algorithms
 - Sleep diaries

35



37

The GGIR pipeline

Part 1	Part 2	Part 3	Part 4	Part 5
- Reads file - Data quality - Extracts features	- Data quality - Descriptive analysis - Reports (csv files)	Heuristic sleep detection	- Sleep analysis - Reports (csv files)	- Full-day time-use analysis - Reports (csv files)

38

The GGIR()

Sleep analysis

GGIR()
[...]
Sleep analysis
HASIBalgo = "vanHees2015",
[...]

GGIR()
[...]
Sleep analysis
HASIBalgo = "ColeKripke1992",
[...]

GGIR()
[...]
Sleep analysis
HASIBalgo = "Sadeh1994",
[...]

GGIR()
[...]
Sleep analysis
HASIBalgo = "Galland2012",
[...]

39

Rest/Wake detection

Algorithms

PLOS ONE

A Novel, Open Access Method to Assess Sleep Duration Using a Wrist-Worn Accelerometer

David E. van Hees^{1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19,20,21,22,23,24,25,26,27,28,29,30,31,32,33,34,35,36,37,38,39,40,41,42,43,44,45,46,47,48,49,50,51,52,53,54,55,56,57,58,59,60,61,62,63,64,65,66,67,68,69,70,71,72,73,74,75,76,77,78,79,80,81,82,83,84,85,86,87,88,89,90,91,92,93,94,95,96,97,98,99,100}

doi: 10.1371/journal.pone.0142533

Fundamental Research

Activity-Based Sleep-Wake Identification: An Empirical Test of Methodological Issues

*Ali Salibi, Katherine M. Sharkey and Mary A. Carskadon
E. P. Bradley Hospital-Brown University, East Providence, Rhode Island, U.S.A.

doi: 10.1093/sleep/17.3.201

Journal of Sleep Medicine

Algorithms for using an activity-based accelerometer for identification of infant sleep-wake states during nap studies

Bethany C. Gault^{1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19,20,21,22,23,24,25,26,27,28,29,30,31,32,33,34,35,36,37,38,39,40,41,42,43,44,45,46,47,48,49,50,51,52,53,54,55,56,57,58,59,60,61,62,63,64,65,66,67,68,69,70,71,72,73,74,75,76,77,78,79,80,81,82,83,84,85,86,87,88,89,90,91,92,93,94,95,96,97,98,99,100}

doi: 10.1093/sleep/15.5.481


Technical Note

Automatic Sleep/Wake Identification From Wrist Activity

*Roger J. Cole, *David F. Kripke, *William Gron, *Daniel J. Mullaney and *J. Christian Gillin
*Department of Psychiatry, University of California, San Diego, and
San Diego Veterans Affairs Medical Center, San Diego, California 92161, U.S.A., and
*Cincinnati Children's Hospital Medical Center, Cincinnati, Ohio 45229, U.S.A.

doi: 10.1093/sleep/15.5.481

40



41

Rest/Wake detection


"vanhees2015"

- Interpretable as lack of posture change and lack of movement, regardless of agreement with neurological sleep
- Angle is a more visual concept than magnitude of acceleration

Angle of z-axis (independent of attachment orientation across bands)

DOI:10.1371/journal.pone.0142533

42



43


Rest/Wake detection

"vanhees2015"


$p(\text{Sleep}) = \max(\text{diff}(\text{angle } z))$

$\leq 5^\circ$: sleep


44




45




46




47



48



49



50

The GGIR()

Sleep analysis

```
GGIR(
  [,]
  # Acceleration metrics
  do.zcx = TRUE, do.zcy = TRUE, do.zcz = TRUE,
  do.brondcounts = TRUE,
  # Sleep analysis
  HASIB.algo = "Sadeh1994",
  Sadeh.axis = "Y",
  [,])
```

 Accelting

53

?

What is the best definition of a sustained inactivity bout in GGIR? *(Single choice)*

- ☐ A period of 5 minutes in which the participant is asleep
- ☐ A period of time in which we observe very low movement or posture change in the accelerometer signal and, therefore, might indicate sleep
- ☐ It is a sedentary bout with especially harmful consequences for health

54

?

What is the best definition of a sustained inactivity bout in GGIR? *(Single choice)*

- ☐ A period of 5 minutes in which the participant is asleep
- ☒ A period of time in which we observe very low movement or posture change in the accelerometer signal and, therefore, might indicate sleep
- ☐ It is a sedentary bout with especially harmful consequences for health

55

?

What algorithms can you use in GGIR to detect rest/wake status? *(Multiple choice)*

- ☐ van Hees 2015, based on the Z-angle
- ☐ Sadeh 1994, based on zero-crossing counts
- ☐ Cole-Kripke 1992, based on zero-crossing counts
- ☐ Galland 2012, based on zero-crossing or ActiGraph counts

56

?

What algorithms can you use in GGIR to detect rest/wake status? *(Multiple choice)*

- ☒ van Hees 2015, based on the Z-angle
- ☒ Sadeh 1994, based on zero-crossing counts
- ☒ Cole-Kripke 1992, based on zero-crossing counts
- ☒ Galland 2012, based on zero-crossing or ActiGraph counts

57


Sleep Period Time



www.accelting.com

58

The GGIR pipeline




59

Accelerometer-based sleep assessment


Free living

- Challenge 2: separate daytime and nighttime (Sleep Period Time)
- Automatic algorithms
- Sleep diaries




61

“Guiders” to guide SPT window detection



62

“Guiders” to guide SPT window detection




63

Guiders

Guider	Definition	Relevant arguments to use it
setwindow		
LSH-12		
HDCZA		
HorAngle		
sleep log		


64

(Fixed) Set window




[van Hees et al. Scientific Reports 2018, doi:10.1038/s41598-018-31266-z]


65




66




67




68



71




72




73

HorAngle algorithm




74

HorAngle algorithm




75

HorAngle algorithm



76

HorAngle algorithm



77

The GGIR()

Sleep analysis

```

GGIR(
  [...]
  # Sleep analysis
  def.noc.sleep = c(23, 7),
  [...])
  
```

→ Guider: setwindows

```

GGIR(
  [...]
  # Sleep analysis
  def.noc.sleep = c(0,
  [...])
  
```

→ Guider: L5+/-12

```

GGIR(
  [...]
  # Sleep analysis
  def.noc.sleep = 1,
  HASPT.algo = "HDCZA",
  [...])
  
```

→ Guider: in-built algorithms



78

The GGIR()

Sleep analysis guided by sleeplog

```

GGIR(
  [...]
  # Sleep analysis
  def.noc.sleep = 1,
  HASPT.algo = "HDCZA",
  loglocation = "C:/mysudy/mysleeplog.csv",
  cold = 1, colim = 1, nrights = 7,
  sleeplosep = "",
  sleepwindowType = "TimeInBed",
  [...])
  
```



79

NOTE: If night is not available, in-built algorithms will be used (`def.noc.sleep` & `HASPT.algo`)

Basic sleeplog

`nnights = 4`

`colid = 1` `coln1 = 2`

ID	OnsetL.n1	Wakeup.n1	OnsetL.n2	Wakeup.n2	OnsetL.n3	Wakeup.n3	OnsetL.n4	Wakeup.n4
01	23:00:00	07:00:00	23:45:00	08:20:00	23:15:00	08:00:00	00:30:00	
02	22:30:00	07:30:00	22:35:00	07:00:00	23:45:00	09:05:00	23:44:00	09:00:00
03	23:45:00	07:10:00	00:02:00	08:30:00	22:50:00	07:25:00	23:00:00	07:38:00
04	00:10:00	09:00:00				07:30:00	00:25:00	09:10:00

As documented in: <https://cranr-project.org/web/packages/GGIR/vignettes/GGIR.html>

80

Advanced sleeplog

`nnights = ..`

`colid = 1` `coln1 = 2`

ID	D1_date	D1_wakeup	D1_inbed	D1_nap_start	D1_nap_end	D1_nonwear_off	D1_nonwear_on	D2_date
01	2022-06-24	07:00:00	23:15:00	15:00:00	15:45:00	13:35:00	14:10:00	2022-06-25
02	2022-06-26	07:30:00	23:45:00			09:05:00	10:30:00	2022-06-27
03	2022-07-24	07:10:00	22:50:00	13:02:00	13:30:00	18:00:00	19:10:00	2022-07-25
04	2022-06-14	09:00:00	00:50:00			20:30:00	21:00:00	2022-06-15

As documented in: <https://cranr-project.org/web/packages/GGIR/vignettes/GGIR.html>

82


Advanced sleeplog

- Date columns → "date"
- Wakeup columns → "wakeup"
- Sleep onset columns → "onset", "inbed", "tobed", "lightsout"
- Napping columns → "nap"
- Nonwear columns → "nonwear"


ID	D1_date	D1_wakeup	D1_inbed	D1_nap_start	D1_nap_end	D1_nonwear_off	D1_nonwear_on	D2_date
01	2022-06-24	07:00:00	23:15:00	15:00:00	15:45:00	13:35:00	14:10:00	2022-06-25

As documented in: <https://cranr-project.org/web/packages/GGIR/vignettes/GGIR.html>


83



84



85



86

Guiders

Summary of guiders

Guider	Definition	Relevant arguments to use it
setwindow	Fixed windows for all nights/participants	
L5+/-12	Midpoint of L5 and surrounding 12h window	
HDCZA	HDCZA algorithm (van Hees 2018)	
HorAngle	HorAngle algorithm intended to detect lying posture	
sleep log	Reported sleep diaries (basic or advanced)	

88

Guiders

Summary of guiders

Guider	Definition	Relevant arguments to use it
setwindow	Fixed windows for all nights/participants	def.noc.sleep = c(23, 7)
L5+/-12	Midpoint of L5 and surrounding 12h window	def.noc.sleep = c()
HDCZA	HDCZA algorithm (van Hees 2018)	def.noc.sleep = 1 HASPTalgo = "HDCZA"
HorAngle	HorAngle algorithm intended to detect lying posture	def.noc.sleep = 1 HASPTalgo = "HorAngle" longitudinal_axis = 2 sensor.location = "hip"
sleep log	Reported sleep diaries (basic or advanced)	log.location = "C:/mystudy/sleeplog.csv" colid = 1, colint = 2, nnights = 7, sleeplogsep = ":", sleepwindowType = "TimeInBed"

89

"Unusual" sleepers

- More than one sleep period time per day?
- Daysleeper?
 - If guider-defined wake-up > 12pm → re-do sleep analysis from 6pm-to-6pm
 - Classified as daysleeper in reports
 - Intended to adapt the algorithm to night workers

90

Other input arguments

```
GGIR(
  [...]
  # Data cleaning
  do.report = 4,
  includenightcrit = 16,
  [...])
```


16 hours available from noon-to-noon or from 6pm-to-6pm

91

Other forms over evaluation...



<https://www.nature.com/articles/s41467-019-09576-1>



<https://www.nature.com/articles/s41467-018-08259-7>

92



What is the sleep period time? (Single choice)

- ☐ The main (longest) episode of rest in the day that starts and ends with sleep
- ☐ The total time in bed
- ☐ The sum of the minutes classified as sleep every day
- ☐ The total sleep time

95



What is the sleep period time? *(Single choice)*

- ☒ The main (longest) episode of rest in the day that starts and ends with sleep
- ☐ The total time in bed
- ☐ The sum of the minutes classified as sleep every day
- ☐ The total sleep time

96



How to instruct GGIR to use the HorAngle algorithm? *(Single choice)*

- ☐ SPT.algo = "HorAngle"
- ☐ HASPT.algo = 2
- ☐ SPT.algo = 2
- ☐ HASPT.algo = "HorAngle"

97



How to instruct GGIR to use the HorAngle algorithm? *(Single choice)*

- ☐ SPT.algo = "HorAngle"
- ☐ HASPT.algo = 2
- ☐ SPT.algo = 2
- ☒ HASPT.algo = "HorAngle"


98

Part 3 & 4 Output

www.accelting.com

102

The GGIR pipeline



103

Part 3 & 4 Output

- results
 - QC
 - part4_nightsummary_sleep_full.csv
 - part4_summary_sleep_full.pdf
 - part4_nightsummary_sleep_cleaned.csv
 - part4_summary_sleep_cleaned.csv
 - visualisation_sleep.pdf

Accelting

104

Part 3 & 4

Output

- results
 - QC
 - part4_nightsummary_sleep_full.csv
 - part4_summary_sleep_full.pdf
 - part4_nightsummary_sleep_cleaned.csv
 - part4_summary_sleep_cleaned.csv
 - **visualisation_sleep.pdf**

 **Accelting**
Accelerating sleep research & sleep medicine


105

Sleep visualizations

do.visual = TRUE

est files > output > output_GTEX > results

Page 1



106

Part 3 & 4

Output

- results
 - QC
 - part4_nightsummary_sleep_full.csv
 - part4_summary_sleep_full.pdf
 - part4_nightsummary_sleep_cleaned.csv
 - part4_summary_sleep_cleaned.csv
 - **visualisation_sleep.pdf**

 **Accelting**
Accelerating sleep research & sleep medicine

108

Output from Part 4

Night-level features (full report)

ID	Onset_n1	Wakeup_n1	Onset_n2	Wakeup_n2	Onset_n3	
01	23:00:00	07:00:00			23:15:00	-

part4_nightsummary_sleep_full.csv

ID	night	cleaningcode	sleeplog_used	acc_available	guider
1	1	0	1	1	sleeplog
1	2	1	0	1	HDCZA
1	3	0	1	1	sleeplog
1	4	0	1	1	sleeplog
1	5	0	1	1	sleeplog
1	6	0	1	1	sleeplog
1	7	2	1	0	sleeplog

Cleaning code	Meaning
0	No problem
1	Sleep log not available (HDCZA used)
2	Not enough valid data
3	No accelerometer data available
4	No nights to be analyzed
5	Guider-defined SPT
6	SPT could not be defined either by sleeplog or algorithms

109

Output from Part 4

Night-level features (full report)

do.report = 4

part4_nightsummary_sleep_full.csv

ID	night	sleeponset	wakeup	SptDuration	SleepDurationSpt	WASO	sleepefficiency	sleeplatency
1	1	27.201	32	4.799	3.696	1.103	0.864	0.215
1	2	26.11	31.936	5.826	4.888	0.938	0.868	0.282
1	3	26.5	32	5.5	4.997	0.503	0.834	0.004
1	4	25	32.475	7.475	7.115	0.36	0.799	0.512
1	5	26.613	32	5.588	5.071	0.517	0.798	0.428
1	6	27.131	32	4.869	3.526	1.343	0.929	0.003
1	7	25.074	31.761	6.688	6	0.688	0.856	0.261

111

Output from Part 4

GGIR [.]
Data cleaning
includenightcrit = 16,
excludelastpart4 = FALSE,
excludelastpart4 = FALSE,
[.]

Night-level features (clean report)

do.report = 4

part4_nightsummary_sleep_cleaned.csv

ID	night	sleeponset	wakeup	SptDuration	SleepDurationSpt	WASO	sleepefficiency	sleeplatency
1	1	27.201	32	4.799	3.696	1.103	0.864	0.215
1	2	26.11	31.936	5.826	4.888	0.938	0.868	0.282
1	3	26.5	32	5.5	4.997	0.503	0.834	0.004
1	4	25	32.475	7.475	7.115	0.36	0.799	0.512
1	5	26.613	32	5.588	5.071	0.517	0.798	0.428
1	6	27.131	32	4.869	3.526	1.343	0.929	0.003
1	7	25.074	31.761	6.688	6	0.688	0.856	0.261

112

Output from Part 4

Night-level features (clean report)

do:report = 4

```
GGIR(
  [...]
  # Data cleaning
  includenightcrit = 16,
  excludelastpart4 = TRUE,
  excludelastpart4 = FALSE,
  [...])
```

part4_nightsummary_sleep_cleaned.csv

ID	night	sleeponset	wakeup	SptDuration	SleepDurationInSpt	WASO	sleepefficiency	sleeplatency
1	1	27.201	32	4.799	3.696	1.103	0.804	0.215
1	3	26.5	32	5.5	4.997	0.503	0.834	0.004
1	4	25	32.475	7.475	7.115	0.36	0.799	0.512
1	5	26.413	32	5.588	5.071	0.517	0.798	0.428
1	6	27.131	32	4.869	3.526	1.343	0.929	0.003

113

Output from Part 4

Night-level features (clean report)

do:report = 4

```
GGIR(
  [...]
  # Data cleaning
  includenightcrit = 16,
  excludelastpart4 = TRUE,
  excludelastpart4 = FALSE,
  [...])
```

part4_nightsummary_sleep_cleaned.csv

ID	night	sleeponset	wakeup	SptDuration	SleepDurationInSpt	WASO	sleepefficiency	sleeplatency
1	1	27.201	32	4.799	3.696	1.103	0.804	0.215
1	3	26.5	32	5.5	4.997	0.503	0.834	0.004
1	4	25	32.475	7.475	7.115	0.36	0.799	0.512
1	5	26.413	32	5.588	5.071	0.517	0.798	0.428
1	6	27.131	32	4.869	3.526	1.343	0.929	0.003

114

Part 3 & 4

Output

- meta
 - sleep.qc
 - graphperday_id_01.pdf
 - ...
- results
 - QC
 - part4_nightsummary_sleep_full.csv
 - part4_summary_sleep_full.pdf
 - part4_nightsummary_sleep_cleaned.csv
 - part4_summary_sleep_cleaned.csv
 - visualisation_sleep.pdf



115

Output from Part 4

Person-level features

do:report = 4

part4_summary_sleep_cleaned.csv

ID	sleeplog_used	n_nights_acc	n_nights_sleeplog	n_WE_nights_complete	n_WD_nights_complete	n_WEnights_daysleeper	n_WDnights_daysleeper
1	1	6	6	2	4	0	0

All days

ID	SptDuration_AD-TSAS_mn	SptDuration_AD-TSAS_sd	SleepDurationInSpt_AD-TSAS_mn	SleepDurationInSpt_AD-TSAS_sd	WASO_AD-TSAS_mn	WASO_AD-TSAS_sd
1	5.372	1.142	4.685	1.375	0.687	0.383

Weekdays

ID	SptDuration_WD-TSAS_mn	SptDuration_WD-TSAS_sd	SleepDurationInSpt_WD-TSAS_mn	SleepDurationInSpt_WD-TSAS_sd	WASO_WD-TSAS_mn	WASO_WD-TSAS_sd
1	5.64	1.277	5.045	1.404	0.594	0.34

Weekend days

ID	SptDuration_WE-TSAS_mn	SptDuration_WE-TSAS_sd	SleepDurationInSpt_WE-TSAS_mn	SleepDurationInSpt_WE-TSAS_sd	WASO_WE-TSAS_mn	WASO_WE-TSAS_sd
1	4.838	0.877	3.966	1.409	0.872	0.532

116

Output from Part 4

Complete list of variables

Sleeponset, sleepset.ts	duration_sib_wakinghours	calendar.date
Wakeup, wakeup.ts	number_sib_sleepperiod	filename
SptDuration	number_of_awakenings	cleaningcode
Sleepparam	number_sib_wakinghours	sleeplog_used
guider_inbedStart, guider_inbedStart.ts	duration_sib_wakinghours_atleast15min	acc_available
guider_inbedEnd, guider_inbedEnd.ts	sleeplatency	guider
guider_inbedDuration	sleepefficiency	SleepRegularityIndex
fraction_night_invalid	page	SnFractionValid
SleepDurationInSpt	daysleeper	longitudinal_axis
WASO	weekday	nonwear_perc_spt

117



What is the column cleaningcode? (Single choice)

- ☐ An indicator of the nights that should be removed from datasets
- ☐ An indicator of the time that the device was not worn during nights
- ☐ An indicator of the quality of the sleep period time detection
- ☐ An estimate of the fraction of time in which the device recorded valid data throughout the measurement

118



What is the column cleaningcode? *(Single choice)*

- ☐ An indicator of the nights that should be removed from datasets
- ☐ An indicator of the time that the device was not worn during nights
- ☒ An indicator of the quality of the sleep period time detection
- ☐ An estimate of the fraction of time in which the device recorded valid data throughout the measurement

119

Assignment 3 (build on assignment 2 script)

1. Run GGIR parts 3 and 4 using:
 - Detect SIBs with the "vanHees2015" algorithm
 - Detect SPT with the "HDCZA" algorithm
2. Look up the output and visualizations

120

Day Evaluation

121

Thank you!



122